

ELITE CURVED OUTDOOR STAIRLIFT



BRUNO ELITE CURVED OUTDOOR STAIRLIFT

WEATHER-TESTED RELIABILITY

Engineered to withstand harsh weather conditions, Bruno's bespoke Elite curved outdoor stairlift offers year-round performance, heavy-duty capacity and durable engineering.

Lifting Capacity
30 st (190 kg)

Temperature Tested
-18°C to 52°C

Learn More
[Bruno.com/UK](https://bruno.com/UK)





ADJUSTABLE COMFORT

Change seat and footrest height and space between armrests to accommodate different body sizes.

WEATHER-READY MATERIALS

Secure all-weather protection with marine-grade vinyl, lasting outdoor paint and a lightweight, waterproof cover.

BESPOKE RAIL

Get the exceptional fit of a handcrafted, vertical rail that maximizes open space. The covered gear rack helps reduce exposure to grease and gives sleek appearance.





SPECIALIZED OPTION

Larger Footrest

Increase footrest size from
254 mm x 229 mm to
305 mm x 279 mm



SAFE, SIMPLE OPERATION

Glide along the stairs by pushing an armrest control paddle.
Rotate seat up to 90° at top of stairs for a safe exit. Seat belt and
obstruction sensors standard.

PEACE-OF-MIND DEPENDABILITY

Gain the confidence of performance testing from -18°C to 52°C,
and battery-powered operation that ensures operation even during
power outage.

5-YEAR LIMITED WARRANTY

Superior engineering backed by a comprehensive warranty.

Learn More

[Bruno.com/UK](https://bruno.com/UK)

Model:
CRE-2110E



TECH SPECS

Lift Capacity:

30 st (190 kg)

Seat Size:

457 mm x 381 mm

Seat Heights*:

552 mm, 577 mm, 603 mm, 629 mm

Seat To Footrest*:

Minimum 480 mm, Maximum 559 mm

Standard Footrest Size:

254 mm x 229 mm

Larger Footrest Option:

305 mm x 279 mm

Curved Rail:

Installs within 190 mm of the wall

*Dimensions are approximate and depending on the slope, riser height and custom seat adjustments made for the user.

ABOUT BRUNO

Veteran-Founded | Family-Owned

Bruno is a premier manufacturer of straight and curved stairlifts. Bruno's award-winning indoor and outdoor stairlifts enable people to access all levels of their homes in comfort. Bruno chairlifts are tested to American, Canadian and European standards to ensure customers travel safely and securely on their stairways.



© Bruno Independent Living Aids, Inc. 2022
2022-08-04-elite-outdoor-curved-uk.indd

For more information please contact:
tvmc mobility on 01242 663800
enquiries@tvmcltd.co.uk