



BRUNO ELITE STRAIGHT OUTDOOR STAIRLIFT

GLIDE UP STAIRS EFFORTLESSLY

Go up and down your outdoor stairs with ease with the heavy-capacity, high-performance Bruno Elite. Ready to use in every season, the Bruno Elite straight outdoor stairlift combines comfort, simplicity and durability.

Lifting Capacity
30 st (190 kg)

Temperature Tested
-18°C to 52°C

Learn More
[Bruno.com/UK](https://bruno.com/UK)

ADJUSTABLE COMFORT

Change seat and footrest height and space between armrests to best match different body types.

WEATHER-RESISTANT MATERIALS

Secure all-weather protection with marine-grade vinyl seat, lasting outdoor painted surfaces and a lightweight, waterproof cover.

COMPACT DESIGN

Gain maximum open space on the stairs with Bruno's vertical rail. Fold arms, seat and footrest when not in use to gain even more space.





SAFE, SIMPLE OPERATION

Glide along the stairs by pushing an armrest control paddle. Swivel seat up to 90° at top of stairs for a safe exit. Seat belt and obstruction sensors standard.

PEACE-OF-MIND DEPENDABILITY

Get the confidence of performance testing from -18°C-52°C and battery-powered operation that ensures operation even during power outage. Backed by 5-year limited warranty.

Learn More
[Bruno.com/UK](https://bruno.com/UK)

Model:
SRE-2010E



TECH SPECS

Lift Capacity:
30 st (190 kg)

Seat Size:
457 mm x 381 mm

Seat Heights*:
552 mm, 577 mm, 603 mm, 629 mm

Seat To Footrest*:
Minimum 480 mm, Maximum 559 mm

Footrest Size:
254 mm x 229 mm

Vertical Rail:
Installs within 152 mm of the wall

*Dimensions are approximate and depending on the slope, riser height and custom seat adjustments made for the user.

ABOUT BRUNO

Veteran-Founded | Family-Owned

Bruno is a premier manufacturer of straight and curved stairlifts. Bruno's award-winning indoor and outdoor stair lifts enable people to access all levels of their homes in comfort. Bruno stairlifts are tested to American, Canadian and European standards to ensure customers travel safely and securely on their stairways.



© Bruno Independent Living Aids, Inc. 2023
2023-03-27-elite-outdoor-straight-uk.indd

For more information please contact:
tvmc mobility on 01242 663800
enquiries@tvmcltd.co.uk