



BRUNO ELITE STRAIGHT STAIRLIFT

RELIABLE PERFORMANCE

Gain high performance, heavy lift capacity and solid dependability with a Bruno Elite straight stairlift. Glide up and down stairs with a smooth, quiet ride and customisable comfort.

Lifting Capacity
30 st (190 kg)

High Performance
Adjustable Comfort

Learn More
[Bruno.com/UK](https://bruno.com/UK)



MULTI-BODY FIT

Standard Adjustability

Change seat and footrest height and space between armrests to match body type.

Larger Option

Standard generous-sized seat. Optional extra-wide seat.

SLEEK RAIL DESIGN

Slim Rail

Vertical rail installs close to wall to leave maximum space on steps.

Covered Gear Rack

Reduces exposure to grease, gives sleek appearance and avoids pinch-points for children and pets.





POWERFUL OPTIONS

Power Footrest

Fold the seat, and the footrest will automatically fold. Ideal for those with back problems or inability to lean over.

Power Swivel

Push an armrest control to swivel seat to exit. Good for people without dexterity to lift seat lever or ability to guide swivel action with feet.

Power or Manual Folding Rail

Fold the bottom rail section to avoid creating an obstruction when nearby doorway or narrow hall.

Wider Seat

Increase size of Elite seat from 457 mm x 381 mm to 572 mm x 419 mm.

SAFE, SIMPLE OPERATION

Easy Operation

Armrest control paddle powers stairlift up and down. Two wireless remotes also included.

Swivel Exit

Seat swivels up to 90° and locks in place for safe exits away from stairs. Seat belt and obstruction sensors further safety.

PEACE-OF-MIND DEPENDABILITY

Continuous Charge

Battery-powered performance even during power outage. Stairlift continuously charges anywhere it stops on rail.

Lifetime Limited Warranty

Superior engineering backed by the industry's most comprehensive warranty.

Learn More

[Bruno.com/UK](https://bruno.com/UK)

Model:
SRE-2010



TECH SPECS

Lift Capacity:
30 st (190 kg)

Standard Seat Size:
457 mm x 381 mm

Larger Seat Option:
572 mm x 419 mm

Seat Heights*:
552 mm, 577 mm, 603 mm, 629 mm

Seat To Footrest*:
Minimum 480 mm, Maximum 559 mm

Standard Footrest Size:
254 mm x 229 mm

Larger Footrest Option:
305 mm x 279 mm

Vertical Rail:
Installs within 152 mm of the wall

Complies with the new EN8140 Regulations.

**Dimensions are approximate and depending on the slope, riser height and custom seat adjustments made for the user. Please note that the addition of a powered swivel seat increases the seat height measurement by 50 mm.*

ABOUT BRUNO

Veteran-Founded | Family-Owned

Bruno is a premier manufacturer of straight and curved stairlifts. Bruno's award-winning indoor and outdoor stairlifts enable people to access all levels of their homes in comfort. Bruno stairlifts are tested to American, Canadian and European standards to ensure customers travel safely and securely on their stairways.



© Bruno Independent Living Aids, Inc. 2023
2023-22-03-elite-straight-uk.indd

For more information please contact:
tvmc mobility on 01242 663800
enquiries@tvmcltd.co.uk